



**Trio Elisse® Sting Free
Skin Barrier**

**say hello
to healthier
skin**

trio Elisse®



do regular bag changes leave your skin constantly red and sore?

Trio Elisse® builds a protective barrier between your stoma bag and your skin.



Prevents skin damage

by stopping output coming in contact with healthy or damaged skin



Protects 'at risk' peristomal skin

from the strain of regular removals by building a barrier between skin and adhesive flange



Allows skin to breathe

naturally during wear leading to better skin health



Soothes irritated skin

Our unique silicone formula calms red and itchy peristomal skin



360° directional spray

sprays at any angle



Dries in seconds

Fast drying, long-lasting protection

Trio Elisse® Silicone Technology

Trio Elisse® Sting Free Skin Barrier has a hypoallergenic 100% silicone formula that means no nasty chemicals and builds a protective barrier between the adhesive flange and your skin. The gentle, directional spray action sprays at any angle without wastage.

how to use Trio Elisse®

Trio Elisse® uses bag in can technology – a recyclable and sustainable packaging system. Free from pollutants used in traditional aerosols, this packaging approach ensures 100% of the product is used.



1. After removing your stoma bag ensure your skin is clean and dry.
2. Hold the can of Trio Elisse® approximately 10cm (4-5 inches) away from skin and spray around your stoma where the adhesive flange will go.
3. Trio Elisse® can be used to build up layers of protection if required.
4. Allow the skin to dry, this should only take a few seconds and then place your stoma bag onto your skin and secure.

Also available in wipes

Handy size, individually wrapped for when you are travelling or out and about.



Trio Elisse® is available on prescription and stocked by all dispensing suppliers and ostomy delivery companies.

We hope you've enjoyed using Trio Elisse®. Please don't hesitate to get in touch with a Trio advisor should you like any further information or need additional advice.

Trio Elite® spray


Code	Pack size
TR204	50ml spray

Trio Elite® wipes

Code	Pack size
TR203	30 wipes



 trioostomycare.com

 0800 531 6280

 careline@triohealthcare.co.uk

