



trio ostomy care™

trio ostomy care™

TRIO HEALTHCARE LTD . MISSENDEN ABBEY . LONDON ROAD  
GREAT MISSENDEN . BUCKINGHAMSHIRE . HP16 0BD

T: 0800 531 6280 E: [mail@triohealthcare.co.uk](mailto:mail@triohealthcare.co.uk)

[TrioOstomyCare.com](http://TrioOstomyCare.com)

---

trusted results in ostomy

Silken®  
silicone stoma gel

Feel the difference, experience the benefit  
of our advanced silicone technology...

trusted results in ostomy



## Silken® – silicone stoma gel

– one piece removal, leaves no residue\*

**Have you ever wished that you could remove your hydrocolloid paste or strip in one piece leaving no residue? We thought you should be able to, so we developed Trio Silken®.**

Trio Silken® is a silicone gel that is used to fill skin folds, creases and uneven areas around the peristomal skin to provide a better surface for the application of an ostomy appliance. Unlike hydrocolloid pastes or strips, Trio Silken® will not absorb moisture or body waste and will not break down, leaving you with a difficult and messy substance to remove. Instead, as wear time progresses, Trio Silken® cures to a single piece of silicone adhesive that removes in one piece when fully cured.

Being sting-free on application and allowing your skin to breathe, Trio Silken® remains soft and flexible to ensure comfort. Trio Silken® can be used to create a perfect custom seal around the stoma, particularly useful for irregular-shaped or difficult to manage stomas.

Benefit from the improved quality of life through the avoidance of effluent leakage and consequential skin excoriation and faecal malodour – try Trio Silken® soft silicone stoma gel and be amazed at how easy it is to remove!

- Fills skin folds, creases and uneven areas around the stoma
- Creates a smooth surface for easy pouch adhesion
- Will not sting on application
- More hygienic than traditional hydrocolloid products – does not absorb faeces or urine
- Remains soft and flexible to ensure comfort
- Allows the skin to breathe
- One piece removal – leaves no residue\*
- Atraumatic removal – less skin irritation



Innovative use of silicone technology to produce a soft silicone stoma gel that fills surgical scars, skin folds, creases and uneven surfaces, removing in one piece when fully cured. Allows the skin to breathe and remains soft and flexible to ensure comfort.

\*when fully cured

Order your **FREE** sample at: [TrioOstomyCare.com](http://TrioOstomyCare.com)