



trio ostomy care™

trio ostomy care™

TRIO HEALTHCARE LTD . MISSENDEN ABBEY . LONDON ROAD
GREAT MISSENDEN . BUCKINGHAMSHIRE . HP16 0BD

T: 0800 531 6280 E: mail@triohealthcare.co.uk

TrioOstomyCare.com

trusted results in ostomy

Siltac®
silicone ostomy seal

Feel the difference, experience the benefit
of our advanced silicone technology...

trusted results in ostomy

trio **Siltac**® – silicone ostomy seal
– skin friendly protection

To date, you have had a limited choice in ostomy seals.
All current products comprise of a hydrocolloid formulation.
All current products suffer from the same set of problems.

Hydrocolloid seals or rings absorb moisture (and body waste!). As the hydrocolloid 'swells' the seal around the stoma becomes better. Level of absorption therefore becomes one of the most important attributes of a hydrocolloid seal.

We believe that having an adhesive containing faeces and urine on your skin for lengthy periods of time is not a hygienic solution to protect against leakage. The removal of a hydrocolloid seal or ring is sometimes difficult due to the absorption of moisture (and body waste!) which leads to the breakdown and failure of the product. This results in the need to repeatedly clean the skin to remove residue and ensure correct fitment of your next ostomy appliance.

Trio **Siltac**® is a soft silicone ostomy seal that works completely differently. It does not absorb moisture or body waste; for the first time you can have a clean and hygienic seal. Trio **Siltac**® seals have a natural tendency to return to their original shape when stretched. This means you simply stretch a Trio **Siltac**® seal around your stoma, and shape to form a perfect custom fit. Trio **Siltac**® removes in one piece leaving no residue – your pouch change regime can now be quicker, easier and cleaner.

- Maintains integrity & shape to provide a longer-lasting seal
- More hygienic than traditional hydrocolloid products – does not absorb faeces or urine
- Allows the skin to breathe
- Reduced risk of leakage – extends wear time of your ostomy appliance
- Stretch around your stoma & shape to create a perfect custom fit
- Trio **Siltac**® conforms to the shape of your stoma
- Promotes comfort due to high level of flexibility
- Instant adhesion – better seal formed
- One piece removal – leaves no residue
- Atraumatic removal – less skin irritation



trio ostomy care™

Innovative use of soft silicone technology to produce an ostomy seal that allows the skin to breathe, stretches to shape and lightly hugs the stoma. Maintains integrity & shape to provide a good seal and extend wear time of ostomy appliances.

Order your **FREE** sample at: TrioOstomyCare.com