



GeniiTM

the stoma bag that
sticks, seals and feels
like it's not even there

trio
GeniiTM



Genii™ sticks, seals and feels



Genii™ is the first ever stoma bag in the world that uses our unique Sil2 Breathable Silicone Technology® in the flange. Which means it sticks with a unique silicone adhesive that flows in and around skin cells allowing stronger adhesion and less risk of leaks.

It seals perfectly to protect the skin from harmful bodily waste, bacteria and moisture and, because the silicone flange moves with your patient's body, it feels like it's not even there.

Genii™ also features a lightweight, water-resistant material which comes in a choice of three skin tone colours.

Feels like it's not even there.

Sticks with a unique silicone adhesive

Flows in and around skin cells meaning stronger adhesion and less risk of leaks

Seals perfectly to protect skin

from harmful bodily waste, bacteria and moisture

Feels like it's not even there

Unique silicone flange moves with the body

Water-resistant material

Leaves the bag dry to the touch. Ideal for being active, swimming, showering and bathing

Skin tone colours

'Camouflages' against the skin - available in a choice of light, medium and dark

Soft material cover

Lightweight, quiet material, stays fresh for optimum, long-lasting comfort

Improved airflow




High performance carbon filter manages air in the bag to reduce ballooning, pancaking and odour

Easy viewing

Split design bag cover to help view contents

Easy to empty and clean (drainable only)

easy opening, draining and cleaning



Genii™ clinical study.¹

Our clinical assessment of Genii™ has been a roaring success.

In a 6 week clinical study of Genii™ we were thrilled to see 35% of ostomates noticed an improvement to their skin, the majority within the first week! Many ostomates reported a better quality of life overall due to improved sleep quality, self-image and social confidence.



Genii™ scores 9/10 for stickability...

“

The soft adhesion disc fitted very well over my hernia and felt very comfortable.

”



Genii™ scores 9.2/10 for comfort...

“

The bag is more comfortable, easy to wear, apply and remove. There has also been less irritation to my skin since using the bag.

”



93% said they would recommend Genii™...

“

It took me a little while to get used to it, but once I did, I wouldn't go back. Just give it time to get used to it – it will be worth it.

”

“

Skin has improved, I'm less self-conscious and not aware of pouch.

”

Genii™ case study:

A noticeable difference to skin.

Patient is a 76-year-old male who underwent surgery for rectal cancer in 2019, the cancer was surgically removed and a stoma created. He had a hernia. Peristomal skin was reported as not flat and mildly sore.

The patient was extremely happy with the Genii™ throughout a 6 week trial. The patient rated it 10 out of 10 for how well it stuck to the skin; comfort compared with his old stoma bag; and security during the night, during exercise and during showering. The positive experience also reflected on certain aspects of the ostomate's life including:



Sleep:

'very much' improved his sleep.
He no longer needs a support belt
and felt it was 'nice and secure'



Exercise routine:

improved 'a little' as he no longer
needed emergency supplies as felt so
confident in the security of the bag



Self-image:

the ostomate also communicated that
Genii™ had improved his self-image
'quite a lot'

There was also a noticeable difference in the ostomate's peristomal skin. The initial photograph taken before the trial shows the skin to be red and there are signs of MARSI (Medical Adhesive Related Skin Injury). A change was noted in the first week and he described 'skin much improved and the skin does not bleed anymore.'



Before Genii™



Week 4

people who've tried Genii™ haven't looked back.

We're confident our new Genii™ stoma bags are going to be life changing. Anne, Alfred and Scott were amongst the first to try it.

Genii™ sticks...

“

What amazed me was the new silicone adhesive stuck much more firmly than anything I've used before. I had no leaks or pancaking. The silicone flange makes all the difference.

Alfred

”

Genii™ seals...

“

I found the silicone product actually clung to my skin. There was no peeling of the flange even after days. It gave me great confidence wearing it and going out to exercise.

Scott

”

Genii™ feels...

“

Once it was on, you couldn't feel it and that's massive because you can forget it's there. It worked incredibly well, there were no gaps and no wrinkling. I think it's a real gamechanger.

Anne

”

combine with a range of Genii™ accessories.

Trio also has a range of essential accessories containing Sil2 Breathable Silicone Technology®. They have all been designed to help your patients combat common problems and promote healthy stoma skin.

Genii™ Ostomy Seals

Genii™ Ostomy Seals protect the delicate skin around the stoma from stool and urine, allowing it to heal naturally and remain healthy. Quick and easy to apply, they simply stretch to provide a custom fit, without the need for warming or cutting.



Genii™ Convex Ostomy Seals

Dips, gaps or uneven skin around the stoma can cause leaks. Genii™ Convex Ostomy Seals are ideal for flat and inverted stomas. They fit smoothly to fill any gaps and apply gentle pressure, if needed, around the stoma.



Genii™ Flange Extenders

If your patient has a hernia, does a lot of physical activity or wants peace of mind at night, Genii™ Flange Extenders offer extra security. They're designed to go around the flange to form a soft, flexible, invisible border which conforms to their body shape.

Genii™ Ostomy Seals

Code	Size	Stoma size	Pack size
TR1020	Small	20–35mm	30
TR1028	Medium	35–50mm	30
TR1035	Large	50–80mm	30

Genii™ Convex Ostomy Seals

Code	Size	Stoma size	Pack size
TR1020CX	Small	20–35mm	10
TR1030CX	Large	30–50mm	10

Genii™ Flange Extenders

Code	Size	Stoma size	Pack size
TR1060	—	All	20

ready to start prescribing Genii™?

Light						
Pip Codes	Code	Size	Flange	Bag size	Filter	Pack Size
Closed bags						
416-5304	TR321N	15-60mm	cut to fit	Mini	✓	30
416-5270	TR322N	15-60mm	cut to fit	Midi	✓	30
416-5247	TR323N	15-60mm	cut to fit	Maxi	✓	30
Drainable bags						
416-5395	TR324N	15-60mm	cut to fit	Mini	✓	30
416-5361	TR325N	15-60mm	cut to fit	Midi	✓	30
416-5338	TR336N	15-60mm	cut to fit	Maxi	✓	30

Medium						
Pip Codes	Code	Size	Flange	Bag size	Filter	Pack Size
Closed bags						
416-5296	TR321T	15-60mm	cut to fit	Mini	✓	30
416-5262	TR322T	15-60mm	cut to fit	Midi	✓	30
416-5239	TR323T	15-60mm	cut to fit	Maxi	✓	30
Drainable bags						
416-5387	TR324T	15-60mm	cut to fit	Mini	✓	30
416-5353	TR325T	15-60mm	cut to fit	Midi	✓	30
416-5320	TR336T	15-60mm	cut to fit	Maxi	✓	30

Dark						
Pip Codes	Code	Size	Flange	Bag size	Filter	Pack Size
Closed bags						
416-5288	TR321C	15-60mm	cut to fit	Mini	✓	30
416-5254	TR322C	15-60mm	cut to fit	Midi	✓	30
416-5221	TR323C	15-60mm	cut to fit	Maxi	✓	30
Drainable bags						
416-5379	TR324C	15-60mm	cut to fit	Mini	✓	30
416-5346	TR325C	15-60mm	cut to fit	Midi	✓	30
416-5312	TR336C	15-60mm	cut to fit	Maxi	✓	30

Bags with split cover

All sizes and shades are available on prescription and stocked by all dispensing suppliers and ostomy delivery companies.

Trio has a team of customer care advisors to help patients with anything they need.

Call us on
0800 531 6280 (UK)
Email us at
careline@triohealthcare.co.uk
Visit
trioostomycare.com

 **Genii™**
Genius for skin.



 trioostomycare.com  0800 531 6280  careline@triohealthcare.co.uk



Trio Healthcare Ltd. Registered in England. No: 05776276. Registered office: Units 2A-2B, Sandylands Business Centre, Carleton New Road, Skipton, North Yorkshire, BD23 2AA, United Kingdom. Sil2 Breathable Silicone Technology® and Trio Genii™ are trademarks of Trio Healthcare Ltd. ©2022 Trio Healthcare Ltd.

M.0163